




Airstream® G4

Class II Type A2 Biological Safety Cabinets

*The Leading Energy-efficient, Quiet,
and Ergonomic Biosafety Cabinet in the World*





USB Port and Zero Volt Relay Contact

- USB Port to send operational information to Building Management System (BMS)
- Zero Volt Relay Contact to turn ON/OFF exhaust blower and signal the building alarm



Centurion Touchscreen Controller

- 7" Capacitive Touch Screen, intuitive to use like phones
- Displays key safety information graphically on one screen
- Built-in guide to use the cabinet and respond to situations
- Centered and angled down for easy reach & viewing
- Datalogger to assist diagnostics and send info to BMS
- Optional: 21 CFR Part 11 Compliance



One-Piece Interior Wall

- Double layer side walls with negative pressure
- Easy to reach service fixture and outlets
- Large radius corner for easy cleaning



Unique Stainless Steel (S-Series) and Glass Hybrid (E-Series) Wall

- Large corner radius for easy cleaning
- Easy to reach service fixture and outlets



Removable Paper Catch

- Easy-to-clean
- Optional pre-filter can be fitted



Raised Ergonomic Arm Rest

- Helps prevent grille blocking
- Comfortable working posture



Dished Work Tray

- Contains spillage
- Angled edge for easy cleaning
- V-shaped grille to prevent blocking
- Work tray holder for drain pan cleaning



Available in 3', 4', 5', and 6'



Airstream® NS G4

Class II Type A2 • Biological Safety Cabinets



Airstream®
Class II Type A2 Biosafety Cabinet



Airflow Sensor

- Monitors real-time airflow for safety
- Alerts the user if airflow is insufficient



Energy-efficient DC ECM Blower

- The leading energy efficient Class II Type A2 Biosafety Cabinet in the world with 70% energy savings compared to AC motor
- Stable airflow despite building voltage fluctuations and filter loading
- Standby mode to further reduce power consumption by 60%




ULPA Filter

- 10x Filtration efficiency of HEPA filter
- Creates ISO Class 3 work zone instead of industry-standard ISO Class 5
- Same 10 years filter life and replacement cost as HEPA filters

Dynamic Chamber™

- Blower plenum and side walls are surrounded by negative pressure
- Prevents contaminants from escaping outside
- Steel plenum with Isocide™ coating, instead of flexible bag plenum, to prevent contamination leak



User-Modified Pass-Through / Cable Port

- 3" port with 1/4" hole on the internal silicone membrane
- NSF-approved and surrounded by negative pressure
- Allows cables and tubes to exit with fully closed sash



Angled Sash

- 10° angled front to optimize user comfort, reduce glare and maximize reach into the work area




ISOCIDE™ Powder Coat

- Silver-ion impregnated powder coat
- Inhibits microbial growth to improve safety
- Prevents the plenum from becoming biohazard landfill

| Standards Compliance and Certifications | | | | |
|---|---|---|---|-------------------------------------|
| Performance | Air Quality | Filtration | Electrical Safety | Electromagnetic Compatibility (EMC) |
| NSF / ANSI 49 | ISO 14644.1, Class 3, Worldwide, US Fed Std 209E, Class 1 USA JIS B9920, Class 3, Japan | EN-1822 (H14), Europe IEST-RP-CC001, USA | UL 61010-1 3rd Ed, USA CSA22.2, No.1010-192, Canada IEC61010-1, Worldwide | EN IEC 61326-1 Group 1 / Class A |

Cabinet Filtration System



Dynamic air barrier, where inflow and downflow converge

Side capture zones

- ULPA-filtered air
- Unfiltered / potentially contaminated air
- Room air / Inflow air


- Ambient air is pulled through front grille to create inflow, without going into the work surface. Inflow is joined by half of the downflow, to create front air curtain that is fine-tuned to create a large performance envelope. The combined air stream travels through the back air column towards the blower.

- Approximately $\frac{1}{3}$ of the air in the common plenum is exhausted through the ULPA filter to the room. The remaining $\frac{2}{3}$ of the air is passed through the downflow ULPA filter and into the work area as a vertical laminar flow air to create ISO Class 3 work surface and prevents cross contamination.

- Near the work surface, the downflow splits. About half goes to the front grille, and half goes to the rear grille. A small portion enters the side capture zones to prevent dead air corners (small blue arrows).


The design was optimized to give large performance envelope, that provides operator and product protection at wide Inflow and Downflow variation from the Nominal point.

The Performance Envelope Concept





- Nominal Airflow
- Personnel / Product Protection
- Area of Personnel and Product Protection
- ▲ NO Personnel / Product Protection
- Area of NO Personnel and Product Protection

Engineering Drawing



- | | | |
|--|---|--|
| 1. Centurion 7" Touchscreen Control System | 7. Service Fixture Retrofit Kit Provision | 13. DC-ECM Blower |
| 2. Tempered Glass - Sliding Sash Window | 8. Electrical Outlet (Max. 5 Amp Total) | 14. Downflow ULPA (H14) Filter |
| 3. Stainless Steel Dished Work Tray | 9. Stainless Steel Side Wall | 15. UV Lamp |
| 4. Electrical/Electronic Panel | 10. Stainless Steel Curved Armrest | 16. Cable Port |
| 5. Dimmable LED lamp | 11. Exhaust Collar (Optional) | 17. Caster wheel - Leveling Feet Combo |
| 6. Provision for IV Bar (Optional) | 12. Exhaust ULPA (H14) Filter | |



Airstream® NS G4

Class II Type A2 • Biological Safety Cabinets

TECHNICAL SPECIFICATIONS

| Model | Stainless Steel Side Walls | 8" Sash Opening | AC2-3S9-NS G4 8" 2011429 | AC2-4S9-NS G4 8" 2011403 | AC2-5S9-NS G4 8" 2011431 | AC2-6S9-NS G4 8" 2011433 |
|--|---|---|---|-------------------------------|-------------------------------|-------------------------------|
| | | 10" Sash Opening | AC2-3S9-NS G4 10" 2011926 | AC2-4S9-NS G4 10" 2011927 | AC2-5S9-NS G4 10" 2011928 | AC2-6S9-NS G4 10" 2011929 |
| | Glass Hybrid Side Walls | 8" Sash Opening | AC2-3E9-NS G4 8" 2011421 | AC2-4E9-NS G4 8" 2011423 | AC2-5E9-NS G4 8" 2011425 | AC2-6E9-NS G4 8" 2011427 |
| Nominal Size | | | 3' | 4' | 5' | 6' |
| External Dimensions (W x D x H) | | | 40.7" x 33.5" x 55.1" | 52.8" x 33.5" x 55.1" | 64.8" x 33.5" x 55.1" | 76.8" x 33.5" x 55.1" |
| Gross Internal Dimensions (W x D x H) | | | 36.0" x 24.0" x 26.0" | 48.0" x 24.0" x 26.0" | 60.0" x 24.0" x 26.0" | 72.0" x 24.0" x 26.0" |
| Usable Work Area (W x D) | | | 35.5" x 19.5" | 47.5" x 19.5" | 59.5" x 19.5" | 71.5" x 19.5" |
| Working Opening | | | 8" / 10" | | | |
| Maximum Opening | | | 19.3" | | | |
| Average Inflow Velocity | | | 105 fpm | | | |
| Average Downflow Velocity | | | 60 fpm | | | |
| Airflow Volume | Inflow / Exhaust Without Ducting | 8" Sash Opening | 210 cfm | 280 cfm | 350 cfm | 417 cfm |
| | | 10" Sash Opening | 263 cfm | 350 cfm | 438 cfm | 525 cfm |
| | Downflow | 8" Sash Opening | 333 cfm | 442 cfm | 554 cfm | 663 cfm |
| | | 10" Sash Opening | 333 cfm | 442 cfm | 554 cfm | 663 cfm |
| | Required Exhaust with Optional Thimble Exhaust Collar | 8" Sash Opening | 224 cfm | 312 cfm | 370 cfm | 459 cfm |
| | | 10" Sash Opening | 274 cfm | 382 cfm | 463 cfm | 559 cfm |
| | Static Pressure For Optional Thimble Exhaust Collar | 8" Sash Opening | 0.06-0.10 in H ₂ O | 0.10-0.14 in H ₂ O | 0.10-0.14 in H ₂ O | 0.18-0.22 in H ₂ O |
| | | 10" Sash Opening | 0.08-0.12 in H ₂ O | 0.12-0.16 in H ₂ O | 0.14-0.18 in H ₂ O | 0.20-0.24 in H ₂ O |
| ULPA Filter Typical Efficiency | | | >99.999% at 0.1 to 0.3 micron, ULPA as per IEST-RP-CC001.3 USA >99.999% at MPPS, H14 as per EN 1822 EU | | | |
| Sound Emission per NSF / ANSI 49* | | | 55.8 dBA | 59.5 dBA | 59.8 dBA | 62.8 dBA |
| LED Lamp Light Intensity | | | 117 ft-cd | 122 ft-cd | 113 ft-cd | 114 ft-cd |
| Cabinet Construction | Main Body | | 0.05" / 18 gauge EG Steel with Isocide™ Oven-Baked Epoxy-Polyester Anti-microbial Powder Coating | | | |
| | Work Zone | | 0.06" / 16 gauge, Stainless Steel 304, 4B finish | | | |
| | Side Walls | Stainless Steel (S Series) | 0.06" / 16 gauge, Stainless Steel 304, 4B finish | | | |
| | | Glass Hybrid (E Series) | 0.25" thickness UV-absorbing Tempered Glass | | | |
| Sash Glass Material | | 0.25" thickness UV-absorbing Tempered Glass | | | | |
| Electrical** 110-130 VAC, 50/60 Hz | Nominal Power Consumption (W) | 8" Sash Opening | 132 | 180 | 242 | 333 |
| | | 10" Sash Opening | 148 | 184 | 255 | 338 |
| | Heat Load (BTU / Hr) | 8" Sash Opening | 450 | 614 | 826 | 1136 |
| | | 10" Sash Opening | 505 | 628 | 870 | 1153 |
| | Cabinet Full Load Amps (FLA) | 8" Sash Opening | 15 | 15 | 15 | 15 |
| | | 10" Sash Opening | 15 | 15 | 15 | 15 |
| Net Weight | | | 448 lbs | 540 lbs | 668 lbs | 796 lbs |
| Shipping Weight | | | 509 lbs | 661 lbs | 814 lbs | 977 lbs |
| Shipping Dimensions, Maximum (W x D x H) | | | 45.3" x 33.5" x 76.8" | 57.1" x 33.5" x 76.8" | 69" x 33.5" x 76.8" | 80.7" x 33.5" x 76.8" |
| Shipping Volume | | | 70.6 cu.ft | 88.2 cu.ft | 105.9 cu.ft | 127.1 cu.ft |


Disclaimer: Technical Specifications may be subjected to further changes without prior notice.

*Noise reading in open field condition / anechoic chamber. Noise reading in normal room varies by room size, layout, and background noise, but may reach roughly 3-4 dBA above these values.

**Electrical power consumption is an measurement of new unit with clean filter operated within nominal setpoint. Result per unit may vary.

Options and Accessories

| | | | | | |
|------------------------------|---|-------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|
| Exhaust Collar | | ECO-F-LA2/AC2-3 G4 5171097 | ECO-F-LA2/AC2-4 G4 5171098 | ECO-F-LA2/AC2-5 G4 5171099 | ECO-F-LA2/AC2-6 G4 5171100 |
| Exhaust Damper | | DAMPER-10 5170104 | | | |
| UV Lamp | | UV-15A-L 5170251 | UV-30A-L 5170255 | | |
| IV Bar | | IV-910 5170499 | IV-1215 5170231 | IV-1520 5170500 | IV-1825 5170501 |
| Electrical Outlet | GFCI | EO-GFCI 5170071 | | | |
| Service Fixtures | US SF-Universal-20 mm | SF-2U22 5170504 | | | |
| | Copper Piping for SF | CU-Pipe 5170026 | | | |
| Support Stand | Telescoping Stand with Caster Wheels | STA-3A0 5131340 | STA-4A0 5131341 | STA-5A0 5131427 | STA-6A0 5131389 |
| | Motorized Stand Height with Levelling Feet | SPML-3A2 5131503 | SPML-4A2 5131504 | SPML-5A2 5131505 | SPML-6A2 5131506 |
| | Motorized Stand Height with Levelling Feet and Seismic Bracket | | SPML-4A2-SB 5131403 | SPML-5A2-SB 5131443 | SPML-6A2-SB 5131404 |
| | Motorized Stand Height with Caster Wheels | SPMC-3A2 5130093 | SPMC-4A2 5130047 | SPMC-5A2 5130100 | SPMC-6A2 5131141 |
| | Motorized Stand with Caster Cradle | | SLC-4A2 G4 with 20" piston 5131441 | SLC-5A2 G4 with 20" piston 5131447 | SLC-6A2 G4 with 20" piston 5131440 |
| Pipette Storage Shelf | | 5260327 | | | |
| Arm Rest Padding | | MEWREST 5170127 | | | |
| Foot Rest | | FT-REST 5170073 | | | |
| Laboratory Chair | | ME-LD-AR360 1150006 | | | |
| IQOQ Protocol | | 9010179 | | | |



ECO-F-LA2/AC2-4 G4



UV LAMP



IV Bar with 6 Hooks



EO-GFCI



SF-2U



STA



SPML



SPMC



SLC



Pipette Storage Shelf



MEWREST



FT-REST



Seismic Bracket



ME-LD-AR360



IQOQ

Life Science Tools



Biological Safety Cabinets



Laminar Flow Cabinet



Animal Workstation



CO2 and Tri-gas Incubators

Chemical Research



Ducted Fume Hoods



Ductless Fume Hoods



cGMP and USP Isolators



Powder Containment Workstations


Pharma Compounding

Medical (IVF)

PCR



Multi-room Incubators



Time-Lapse Incubators



IVF Workstations



PCR Workstations

Bioprocessing

Pharmaceutical



Bioreactors



cGMP and Filling Line Isolators



Air Shower and Pass Boxes



Downflow Booths

ESCO LIFSCIENCES GROUP

42 LOCATIONS IN 24 COUNTRIES ALL OVER THE WORLD



- 📍 Global Offices
- Distributors
- Manufacturing Plants
- R&D Centres
- Regional Distribution Centres

Follow us on social media, download our apps,
and scan the QR code for more info.



@EscoLifesciences



@EscoLifesciences



@EscoLifesci



@Esco



@EscoLifesciences



@EscoLifesciences



Esco Lifesciences



Esco Lifesciences

ESCO

LIFSCIENCES GROUP

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA
Tel: +1 215-441-9661 • Fax 484-698-7757
eti.admin@escolifesciences.com

Esco Micro Pte. Ltd. • 19 Changi South Street 1 • Singapore 486779
Tel +65 6542 0833 • mail@escolifesciences.com
www.escolifesciences.com

Esco Lifesciences Group Offices: Bangladesh | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania | Malaysia | Myanmar | Philippines | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam

9010510_AC2 Gen 4_US_vB_032626

Esco can accept no responsibility for possible errors in catalogues, brochures and other printed materials. Esco reserves the right to alter its products and specifications without notice. All trademarks and logotypes in this material are the property of Esco and the respective companies.