



*Airstream® Horizontal Laminar Flow
Model: LHG-4DS-F_*

Airstream®

Horizontal Laminar Flow Cabinets
Ideal for Plant Tissue Culture



ESCO
SCIENTIFIC



Welcome to Esco

Esco's Vision is to provide enabling technologies for scientific discoveries to make human lives healthier and safer.

The Esco Lifesciences Group is committed to deliver innovative solutions for the clinical, life sciences, research, industrial, laboratory, pharmaceutical, and IVF community. With the most extensive product line in the industry, Esco have passed a number of international standards and certifications. Esco represents innovation and forward-thinking designs, that are of the highest standard quality since 1978.

Availability and Accessibility. Esco has headquarters in Singapore, Indonesia, and Philippines, with manufacturing facilities are located in Asia and Europe. Research and Development (R&D) is conducted worldwide spanning the US, Europe and Asia. Sales, services and marketing subsidiaries are located in 42 major markets including US, UK, Japan, China and India. Esco regional distribution centers are located in Singapore, Malaysia, Thailand, Vietnam, Myanmar, Indonesia, Philippines, Bangladesh, Hong Kong, Taiwan, South Korea, China, Japan, India, UAE, Central and South Africa, Denmark, Germany, Italy, Lithuania, Russia, United Kingdom, and USA. Because of our worldwide presence, you can be sure that Esco is within your reach.

High Quality, Reliable, and Dependable. Esco products are of high quality, reliable, and dependable; assuring customers of research accuracy. Cross functional teams from Esco Production, R&D, Quality Assurance, and Senior Management, are regularly assembled to review and implement areas for improvement.

Esco Cares for Your Safety. Esco focuses on providing safety not just for your samples but also for you and the environment.

Esco Cares for Your Comfort. Building ergonomic designs and reducing noise levels of the units ensures comfort for our users.

Esco Cares for the Environment. One in every four of Esco's employees is involved in R&D and a number of them evaluate new components and/or designs to produce energy efficient equipment. Being GREEN is more than just modifying parts used to produce a new energy efficient technology, it is also embodied in the every aspect of the company.

Customer Service and Support. Our service does not stop once purchase has been done. Esco gives on-time customer service and offers end-user seminars, service training, preventive maintenance, and provides educational materials and informative videos.

As Esco takes the opportunity to respond to the world's needs, we aim not only to contribute in the advancement of scientific discoveries but also in making the world a safer, healthier, and better place to live in.

Airstream®

Laminar Flow Cabinets

Products and Application

Laboratory Equipment

Sample Handling and Preparation

- Class I Biological Safety Cabinets
- Class II Biological Safety Cabinets
- Class II Type A2 Biological Safety Cabinets
- Class II Type B1 Biological Safety Cabinets
- Class II Type B2 Biological Safety Cabinets
- Class III Biological Safety Cabinets
- Horizontal Laminar Flow Cabinets
- Vertical Laminar Flow Cabinets
- Laboratory Animal Research Workstations
- Laboratory Centrifuges

Sample Cultivation

- CO₂ Incubators, Direct Heat Air-Jacketed
- CO₂ Incubators with Cooling System
- CO₂ Incubators with High Heat Sterilization
- Laboratory Shakers

Amplification and Detection

- Conventional Thermal Cyclers
- Microplate Shakers
- PCR Cabinets

Sample Storage & Sample Protection Solutions

- Ultra-low Temperature Freezers
- Lab Refrigerators and Freezers
- Sample Database Management Software
- Intelligent Remote Monitoring Application Protocol
- Remote Monitoring, Datalogging, Programming Software
- Wireless Monitoring System

Chemical Research

- Ducted Fume Hoods
- Ductless Fume Hoods
- Filtered Storage Cabinets
- Powder Weighing Balance Enclosure
- Exhaust Blowers
- Fume Hood Airflow Monitor

General Equipment

Laboratory Thermostatic Products

- Forced Convection Laboratory Oven
- Forced Convection Laboratory Incubator
- Natural Convection Laboratory Incubator
- Refrigerated Laboratory Incubator

Medical / IVF Equipment

Controlled Embryo Handling

- Esco Multi-Zone ART Workstation
- Esco Multi-Zone ART Workstation Class II
- AVT Anti-Vibration Table
- Semi-Closed Environment (SCE) IVF

Safe Embryo Culture

- MIRI® Multiroom Incubator
- MIRI® II Multiroom Incubator
- Mini MIRI® Humidified Incubator
- Mini MIRI® Dry Incubator
- CelCulture® CO₂ Incubator

Innovative Time-Lapse Imaging

- MIRI® Time-Lapse Incubator

Accurate Quality Control

- MIRI® GA Gas and Temperature Validation Unit

Unique Consumables

- CultureCoin®

Healthcare

Esco Pharma Products

Airflow Containment

- BioBooth™
- Ceiling Laminar Airflow (CLAF)
- Cytoculture® Cytotoxic Safety Cabinet
- Pharmacon™ Downflow Booth
- Esco Garment Storage Cabinet
- Esco Glassware Hoods
- Laminar Flow Horizontal/Vertical Trolley (LFH/VT)
- Laminar Flow Straddle Units
- Evidence Drying Cabinet

Isolation Containment

- Advanced Processing Platform Isolator (APPI)
- Aseptic Containment Isoaltor (ACTI)
- Blood Cell Labelling Isolator
- Streamline® Closed Restricted Access Barrier System (SLC-RABS)
- Containment Barrier Isolator (CBI)
 - CBI-Unidirectional (CBI-U)
 - CBI-Turbulent (CBI-T)
 - CBI-Class III Biosafety Cabinet (CBI-III)
 - CBI-Convertible Class III/Class I Biosafety Cabinet (CBI-H)
- Isoclean® Healthcare Platform Isolator (HPI)
 - HPI-G3-Without Filter Below Work Zone
 - HPI-G3-With Filter Below Work Zone
 - HPI-Inflatable Seal (HPI-IS)
- General Processing Platform Isolator
 - GPPI-Inflatable Seal (GPPI-IS)
 - GPPI-Static Seal (GPPI-SS)
- Streamline® Compounding Isolator
 - SCI - Isolator Configuration
 - SCI - Class III Biosafety Cabinet (SCI-III)
- Technetium Dispensing Isolator
- Turbulent Flow Aseptic Isolator (TFAI)
- Weighing and Dispensing Containment Isolator (WDCI)

Cross Contamination Facility Integrated Barrier

- BioPass™ Pass Through
- Cleanroom Air Showers
- Dynamic Pass Boxes/ Dynamic Floor Laminar Hatches
- Infinity® Air Shower Pass Box
- Esco Sputum Booth
- Infinity® Pass Boxes
- Infinity® Cleanroom Transfer Hatch
- Soft Capsule® Soft Wall Cleanroom

Ventilation Containment

- Ventilated Balance Enclosure

Esco VacciXcell Products

Bioreactors and Fermenters

- CelXrock™
- CelCradle™
- CelShaker™
- CelCradle™ X
- CelCradle Semi-Automated Harvesting System™ (CCX-SAH)
- BioXcell™
- StirCradle™
- StirCradle™ PRO
- TideXcell™
- TideXcell™ Harvesting System (TXLHS)
- VXL™ Hybrid Bioreactor

Cell Culture Monitoring, Media and Consumables

- Super Plus™
- Plus™ Vero
- Plus™ MDCK
- Plus™ MDCK II
- BioNOC™ II macrocarriers
- GlucCell™ Glucose Monitoring System
- CVD Kit

Filling Line Equipment

- Filling Line Isolators
- cRabs (close restricted access barriers)
- oRabs (open restricted access barriers)

Integrated Solutions

- Cell Processing Isolator
- Cell Processing Center

Esco TaPestle Rx Products

Pharmacy Compounding Solutions

- Compounding Pharmacy Isolators (SCI, HPI, CBI, GPPI)
- Safety Cabinets and Enclosures (CYT, Class II BSC, VBE, LFC)
- Aseptic Filling Systems


Radiopharmacy Equipment Solutions

- Radioisotope Fume Hood
- Lead-lined Biosafety Cabinet
- Technetium Dispensing Isolator
- Blood Cell Labeling Isolator
- GMP-compliant Radioisotope Dispensing Isolator



Simple Switches Controller

- Easy-to-use and robust
- Pressure gauge to indicate filter load



Isocide™ Antimicrobial Coating

- Silver-ion impregnated powder coat
- Inhibits microbial growth to improve safety




Glass Enclosed Work Zone

- Comfortable unobstructed view
- No electrical outlets, service fixtures and IV bar holes, or gaps, to facilitate easy cleaning and prevent harboring contamination



Single Piece Tray

- Easy-to-clean
- Optional hole for glass bead sterilizer



ULPA Filter


- 10x filtration efficiency than of HEPA filter, creates an ISO Class 3 work-zone instead of the industry standard ISO Class 5
- 10x cleaner work-zone than cabinets with HEPA filter



Airstream® Horizontal Laminar Flow
Model: LHG-4DS-F_

Energy-efficient DC ECM Blower

- 70% Energy savings from AC motor
- Stable airflow, despite building voltage fluctuations
- Night Setback to further reduce power consumption by 60%



Designed to be the Ultimate Green and Serene

Plant Tissue Culture is one of the most promising branches of plant science, occupying key positions in plant biotechnology, plant breeding, and plant propagation. Plant tissue culture, by definition, is the application of various techniques for the aseptic culture of plant cells, tissues, organs, and their components under known physical and chemical conditions *in vitro*.


The remarkable success of the applications of plant tissue culture in plant biotechnology and plant science has led it to its significant roles in agriculture, forestry, and horticulture. Several of its applications include large-scale micro-propagation, somatic embryogenesis and synthetic seed technology, genetic engineering of plants by gene transfer technology, and commercial production of plants.

Among the many concerns in plant tissue culture, the sensitivity of its samples to any source of contamination highly influences its success. A single bacterial cell or fungal spore can be fatal to a small plant piece since it can easily contaminate the growth media. It is essential to maintain a sterile environment as much as possible to avoid the contamination of samples, and this can be provided by the use of a laminar flow cabinet.


The Role of Esco in Plant Tissue Culture

Since 1978, Esco has installed tens of thousands of laminar flow cabinets for the global life sciences market, providing reliable protection for samples and work processes of multitude of applications. Esco laminar flow cabinets are the premium selection for the discerning researchers, offering a combination of value, high quality construction, low operating noise levels and a wide product range.

Introducing Airstream® Horizontal Laminar Flow Cabinets, designed specifically for plant tissue culture. It offers a combination of smooth and easily cleanable work zone, easy-to-use features, superior filtration efficiency, and energy-efficient technology as part of Esco's advocacy to produce environment-friendly laboratory equipment. Experience quality protection. Experience innovation.




Green - Very Low Energy Consumption



Serene - One of the World's Quietest Laminar Flow Cabinet

HORIZONTAL LAMINAR FLOW CABINET ENGINEERING DRAWING



1. Pressure Gauge

2. Operating Switches

3. Stainless Steel Diffuser

4. Stainless Steel Work Tray

5. Stainless Steel Tray

6. Support Stand
7. Pre-filter



8. DC-ECM Blower

9. Tempered Glass Ceiling

10. LED lamp

11. ULPA Filter

12. Tempered Glass Side Panel



Guide to Models, Horizontal Laminar Flow Cabinet

| L H G - D S - F | | | | | | | | | | | | | | | |
|------------------------------|---|----------------------|---|---------------------------|---|-----------------------|---|---|---|-------------------------|---|------------------------|---|----------------------------|-------------------|
| 1st Placeholder Product Line | | 2nd Placeholder Flow | | 3rd Placeholder Side Wall | | 4th Placeholder Width | | 5th Placeholder Internal Width | | 6th Placeholder Control | | 7th Placeholder Window | | 8th Placeholder Electrical | |
| Laminar Flow | L | Horizontal | H | Glass | G | 4 feet | 4 | 3 ft. width, Glass ceiling construction | D | Simple Switches | S | Fixed | F | 230 VAC, 50/60 Hz | 8 |
| | | | | | | | | | | | | | | | 115 VAC, 50/60 Hz |


Example Code: LHG-4DS-F8
Specifications: Laminar Flow Cabinet, Horizontal Flow, Glass sided walls, 4 ft. width, Simple Switches controller, Fixed window, and 230 VAC, 50/60 Hz

| TECHNICAL SPECIFICATIONS | | |
|---|------------------------------|--|
| Model | 230 VAC 50/60 Hz | LHG-4DS-F8 (2120725) |
| | 120 VAC 50/60 Hz | LHG-4DS-F9 (2120726) |
| Nominal Size | | 1.2 meter (4') |
| External Dimensions (W x D x H) | | 1300 x 1065 x 1105 mm (51.2" x 41.9" x 43.5") |
| Gross Internal Dimensions (W x D x H) | | 1185 x 625 x 575 mm (52.8" x 24.5 x 22.6") |
| Usable Work Area | | 0.74 m ² (7.96 ft ²) |
| Airflow | | 0.45 m/s (90 fpm) |
| ULPA Filter Typical Efficiency | | >99.999% at 0.1 to 0.3 micron, ULPA as per IEST-RP-CC001.3 USA >99.999% at MPPS, H14 as per EN 1822 EU |
| Sound Emission per IEST-RP-CC002.2 / NSF 49 | | 52 dBA |
| LED Lamp Intensity (lux) | | >1800 lux (>167 foot candles) |
| Cabinet Construction | Main body | 1.2 mm (0.05") /18 gauge EG Steel With Isocide™ Oven-Baked Epoxy-Polyester Powder Coating |
| | Work Zone | Tray: 1.2 mm (0.08") / 18 gauge, SS 304, 4B Finish |
| | Side Walls | 8 mm (0.3") UV Absorbing Tempered Glass |
| Electrical 230 VAC, 50/60 Hz | Cabinet Full Load Amps (FLA) | 6 |
| | Heat Load (BTU / Hr) | 460 |
| | Nominal Power Consumption | 135 W |
| Electrical 120 VAC, 50/60 Hz | Cabinet Full Load Amps (FLA) | 3.5 |
| | Heat Load (BTU / Hr) | 444 |
| | Nominal Power Consumption | 130 W |
| Net Weight * | | 200 Kg (441 lbs) |
| Shipping Weight * | | 250 Kg (551 lbs) |
| Shipping Dimensions, Maximum (W x D x H)* | | 1550 x 1300 x 1300 mm (45.3" x 51.2" x 51.2") |
| Shipping Volume, Maximum * | | 2.3 m ³ (81.2 ft ³) |

*In an open field condition


ESCO LIFESCIENCES GROUP

42 LOCATIONS IN 21 COUNTRIES ALL OVER THE WORLD



-  Global Offices
-  Licensee
-  Distributors
-  Factories
-  R&D Centers
-  Regional Distribution Centers

Join us on Social Media and Download our Apps!



Filtra Check App
The World's First Chemical
Evaluation App



Esco Social Media
and Apps

ESCO

LIFESCIENCES GROUP

Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777
Tel: +65 6542 0833 • Fax: +65 6542 6920 • mail@escolifesciences.com
www.escolifesciences.com

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA
Tel: +1 215-441-9661 • Fax 484-698-7757
eti.admin@escolifesciences.com

Esco Lifesciences Group Offices: Bangladesh | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam