



Ascent™ Storage Cabinet, Model: ASC-4

Ascent™

Storage Cabinet

Filtered Storage Solution





Welcome to Esco Lifesciences

Esco Lifesciences' vision is to provide enabling technologies for scientific discoveries to make human lives healthier and safer.

Esco Lifesciences is committed to delivering innovative solutions for the clinical, life sciences, research, industrial, laboratory, pharmaceutical, and IVF communities. With the most extensive product line in the industry, Esco has passed a number of international standards and certifications. Esco Lifesciences represents innovation and forward-thinking designs, that are of the highest standard quality since 1978.

Availability and Accessibility. Esco Lifesciences has headquarters in Singapore, Indonesia, and Philippines, with manufacturing facilities located in Asia and Europe. Research and Development (R&D) is conducted worldwide spanning the US, Europe and Asia. Sales, services, and marketing subsidiaries are located in 42 major markets including US, UK, Japan, China and India. Esco regional distribution centers are located in Singapore, Malaysia, Thailand, Vietnam, Myanmar, Indonesia, Philippines, Bangladesh, Hong Kong, Taiwan, South Korea, China, Japan, India, UAE, Central and South Africa, Denmark, Germany, Italy, Lithuania, Russia, United Kingdom, and USA. Because of our worldwide presence, you can be sure that Esco is within your reach.

High Quality, Reliable, and Dependable. Esco Lifesciences products are of high quality, reliable, and dependable. Cross-functional teams from Esco Production, R&D, Quality Assurance, and Senior Management, are regularly assembled to review and implement areas for improvement.

Esco Lifesciences Cares for Your Safety. Esco Lifesciences focuses on providing safety not just for your samples, but also for you and the environment.

Esco Lifesciences Cares for Your Comfort. Building ergonomic designs and reducing noise levels of the units ensure comfort for our users.

Esco Lifesciences Cares for the Environment. Esco Lifesciences incorporates the latest proven technologically advanced components available. One in every four of Esco's employees is involved in Research and Development and are evaluating new components or designs for better efficiency. Whenever a new technology is available, Esco Lifesciences redesigns technology into our new products that will use lesser energy.

Customer Service and Support. Our service does not stop once purchase has been done. Esco Lifesciences gives on-time customer service such as service training, preventive maintenance, and re-certification, to respond to your equipment needs. Esco Lifesciences also offers free end-user seminars and provides educational materials and informative videos.

As Esco Lifesciences takes the opportunity to respond to the world's needs, we aim not only to contribute to the advancement of scientific discoveries but also in making the world a safer, healthier, and better place to live in.



Products and Applications

Life Sciences Laboratory Equipment

Sample Preparation

- Class I Biological Safety Cabinets
- Class II Biological Safety Cabinets
- Class II Type A2 Biological Safety Cabinets
- Class II Type B1 Biological Safety Cabinets
- Class II Type B2 Biological Safety Cabinets
- Class III Biological Safety Cabinets
- Horizontal Laminar Flow Cabinets
- Vertical Laminar Flow Cabinets
- Laboratory Animal Research Workstations
- Laboratory Centrifuges

Sample Cultivation

- CO₂ Incubators, Direct Heat Air-Jacketed
- CO₂ Incubators with Cooling System
- CO₂ Incubators with High Heat Sterilization
- Laboratory Shakers

Amplification and Detection

- Conventional Thermal Cyclers
- Microplate Shakers
- PCR Cabinets

Sample Storage & Sample Protection Solutions

- Laboratory Refrigerators and Freezers
- Ultra-low Temperature Freezers
- Remote Monitoring, Datalogging, and Programming Software
- Wireless Monitoring System

Chemical Research

- Ducted Fume Hoods
- Ductless Fume Hoods
- Filtered Storage Cabinets
- Powder Weighing Balance Enclosure
- Exhaust Blowers
- Fume Hood Airflow Monitor

General Equipment

Laboratory Thermostatic Products

- Forced Convection Laboratory Oven
- Forced Convection Laboratory Incubator
- Natural Convection Laboratory Incubator
- Refrigerated Laboratory Incubator

Medical / IVF Equipment

Innovative Time-Lapse Imaging

- MIRI® TL6 and MIRI® TL12

Embryo Culture

- MIRI® Multiroom Incubator
- MIRI® Humidity Multiroom Incubator
- Mini MIRI® Dry and Humidity Incubator
- CelCulture® CO₂ Incubator
- MIRI® II-12 Multiroom Incubator

Sample Handling

- Esco Multi-Zone ART Workstation
- MIRI® AVT
- Versati™ Tabletop Centrifuge
- Airstream® Laminar Flow Bench

Accurate Quality Control

- MIRI® GA Gas and Temperature Validation Unit

Unique Consumables

- CultureCoin®

Traceability Tool

- MIRI® Evidence

Healthcare

Esco Pharma Products

Airflow Containment

- BioBooth®
- Ceiling Laminar Airflow (CLAF)
- Cytoculture® Cytotoxic Safety Cabinet (CYT)
- Pharmacon™ Downflow Booth
- Esco Garment Storage Cabinet
- Esco Glassware Hoods
- Laminar Flow Horizontal/Vertical Trolley (LFH/VT)
- Laminar Flow Straddle Units
- Evidence Drying Cabinet

Isolation Containment

- Advanced Processing Platform Isolator (APPI)
- Aseptic Containment Isolator (ACTI)
- Blood Cell Labelling Isolator
- Streamline® Closed Restricted Access Barrier System (SLC-RABS)
- Containment Barrier Isolator (CBI)
 - CBI-Unidirectional (CBI-U)
 - CBI-Turbulent (CBI-T)
 - CBI-Class III Biosafety Cabinet (CBI-III)
 - CBI-Hybrid (CBI-H)
- Isoclean® Healthcare Platform Isolator (HPI)
 - HPI-G3-Without Filter Below Work Zone
 - HPI-G3-With Filter Below Work Zone
 - HPI-Inflatable Seal (HPI-IS)
 - HPI - Inflatable Seal - BioVap™ (HPI-IS-BVP)
- General Processing Platform Isolator (GPPI)
 - GPPI-Inflatable Seal (GPPI-IS)
 - GPPI-Static Seal (GPPI-SS)
- Streamline® Compounding Isolator (SCI)
- Streamline® Containment Isolator - Class III (SCI-III)
- Technetium Dispensing Isolator
- Turbulent Flow Aseptic Isolator™ (TFAI™)
- Weighing and Dispensing Containment Isolator (WDCI)

Cross Contamination Facility Integrated Barrier

- BioPass™ Pass Through
- Cleanroom Air Showers
- Dynamic Pass Boxes/ Dynamic Floor Laminar Hatches
- Infinity® Air Shower Pass Box
- Esco Sputum Booth
- Infinity® Pass Boxes
- Infinity® Cleanroom Transfer Hatch
- Soft Capsule® Soft Wall Cleanroom

Ventilation Containment

- Ventilated Balance Enclosure

Esco VacciXcell Products

Tide Motion Bioreactors

- CelXrocker™ (CXR)
- MiniTide®
- CelCradle™ (CC)
- CelCradle X® (CCX)
- TideXcell® (TXL)

Stirred Tank Bioreactors

- BioXcell®
- StirCradle™
- StirCradle™ PRO
- VXL™ Hybrid

Harvesting System

- CelShaker™
- CelCradle X® Semi-automated Harvester System (CCX-SAH)
- TideXcell® Harvester System (TXLHS)

Cell Culture Monitoring, Media and Consumables

- Super Plus™
- Plus™ Vero
- Plus™ MDCK
- Plus™ MDCK II
- BioNOC™ II macrocarriers
- GlucCell™ Glucose Monitoring System
- CVD Kit

Filling Line Equipment

- Traditional Filling Line
- cRABS (Closed Restricted Access Barrier System)
- oRABS (Open Restricted Access Barrier System)

Integrated Solutions

- Cell Processing Isolator (CPI)
- Cell Processing Center

Esco TaPestle Rx Products and Services

Pharmacy Automation and Compounding Supply

- Compounding Pharmacy Isolators (SCI, HPI, CBI, GPPI)
- Safety Cabinets and Enclosures (CYT, Class II BSC, VBE, LFC)
- Filling Line Isolators/RABS

Radiopharmacy Equipment


- Blood Cell Labeling Isolator
- Cytoculture® Lead-Shielded Class II Biosafety Cabinet
- Frontier® Radioisotope™ Fume Hood
- GMP-compliant Radioisotope Dispensing Isolator
- Radiopharmacy Hood (Lead Shielded Biological Safety Cabinet)
- Technetium Dispensing Isolator

Ascent™ Storage Cabinet

INTRODUCTION

Fumes from chemical containers often result to the storage cabinet's corrosion and contaminated air in the laboratory. Ascent™ storage cabinet is ideal in providing a safe and convenient storage solution. This cabinet is equipped with Nanocarb™ filters that effectively adsorb chemical fumes to provide operator and environment protection.

EACH ASC COMES STANDARD WITH:



VOC SENSOR



SENTINEL™ CONTROL SYSTEM



PP TRAYS



HEPA FILTER
(FOR B AND D SERIES ONLY)



FEATURES:



Easy Filter Change



Safety Lock to restrict access

Designed for Enhanced Usability and Efficiency

The Ascent™ Storage Cabinet is designed specifically to provide the operator with a high level of usability, comfort and visibility.

- Industrial-grade main body constructed of electro-galvanized steel: with an abrasion resistant white oven-baked powder-coated finish.
- Esco Isocide™ antimicrobial surface on all painted surfaces minimizes surface contamination.
- Multiple filter configuration options to suit customer's needs.
- Permanently lubricated direct drive centrifugal fan(s); energy efficient, extremely low noise and vibration levels due to proprietary construction and mounting technology.
- Acrylic Glass doors with alarms to notify if doors are left opened.
- Large tray capacity to accommodate more chemical bottles with different weights.

Enhanced Filtration System

Esco's Nanocarb™ activated carbon filters are constructed in order to ensure maximum filter efficiency, retention capacity and operator protection. Esco's research scientists and engineers, working in consultation with world-leading authorities on adsorption science, have developed the following set of unique technologies:


- Optimized retention capacity.
- Generously sized filters with more activated carbon by weight retain more chemicals and last longer.
- Proprietary computer modelling software to predict application suitability, filter saturation capacity, and efficiency for single and multiple compounds.
- Quick-change out filter clamping mechanism allows filter replacements to be carried out with minimal tools; even filter clamping (perimeter, not point force) prevents leaks from occurring.
- Diffusion technology to ensure even filter loading and better airflow uniformity.
- Can filter most commonly used solvents specially acids and organic solvents

Sentinel™ Microprocessor Control System

Esco's Sentinel™ microprocessor control system supervises the operation of all cabinet's functions. It is fully configurable and is equipped with enhanced features to promote cabinet usability. It also has a bright, easy to read LCD panel located at the center of the cabinet.

| Guide to Models | | | | | | | |
|-------------------------------|------|----------------|------|-------------------------------|------|--------------------|------|
| Model | Code | External Width | Code | Filter Configuration | Code | Electrical | Code |
| Ascent™ Storage Cabinet | ASC | 27.0" (2ft.) | 2 | Carbon filter only | A | 230 VAC , 50/60 Hz | 8 |
| | | 53.9" (4ft.) | 4 | Carbon Filter + HEPA Filter | B | 115 VAC , 50/60 Hz | 9 |
| | | 80.9" (6ft.) | 6 | Carbon filter + Carbon filter | C | | |
| | | | | HEPA filter only | D | | |
| | | | | HEPA Filter + Carbon Filter | E | | |

Engineering Drawing



1. Sample Port
2. Sentinel™ Control System

3. Electrical panel
4. Front panel

5. Secondary Filter
(For B, C, and E series only)

6. Blower
7. Main Filter


General Specifications, Ascent™ Storage Cabinet (2ft. series)

| Models | 220-240 VAC, 50-60 Hz, 1Ø | ASC-2A8 2040293 | ASC-2B8 2040294 | ASC-2C8 2040295 | ASC-2D8 2040296 | ASC-2E8 2040328 |
|---------------------------------|-------------------------------|--|-----------------------------|-----------------------------|--------------------|-----------------------------|
| | 110-120 VAC, 50-60 Hz, 1Ø | ASC-2A9 2040297 | ASC-2B9 2040298 | ASC-2C9 2040299 | ASC-2D9 2040300 | ASC-2E9 2040330 |
| Nominal Size | | 0.7 m (2') | | | | |
| External Dimensions (W x D x H) | | 685 mm x 700 mm x 2050 mm (27.0" x 27.6" x 80.7") | | | | |
| Door | | Powder Coated electro-galvanized steel door with 6 mm thick Acrylic window and lockable Swing Handle | | | | |
| Main Body | | 1.2 mm 18 gauge electro-galvanized steel with white oven-baked epoxy-polyester ISOCIDE™ antimicrobial powder coated finish | | | | |
| Number of Blower/s | | 1 | | | | |
| Air volume | | 140 cfm | | | | |
| Trays | Number of Trays | Standard: 4 Trays Maximum: 6 Trays | | | | |
| | Tray capacity per bottle size | For 2.5 - 4 L bottles: each tray can fit 9 bottles For 2.0 L bottles: each Tray can fit 12 bottles | | | | |
| | Maximum weight per tray | 40 kg | | | | |
| | Tray dimension | 600 mm x 600 mm (23.6" x 23.6") | | | | |
| | Tray material | Polypropylene | | | | |
| Filtration Elements | Main Filter | Nanocarb™ Filter (type A-H) | Nanocarb™ Filter (type A-H) | Nanocarb™ Filter (type A-H) | HEPA filter | HEPA Filter |
| | Secondary Filter | None | HEPA Filter | Nanocarb™ Filter (type A-H) | None | Nanocarb™ Filter (type A-H) |
| Control System | | Sentinel™ Silver microprocessor control system | | | | |
| VOC Sensor | | ✓ | | | | |
| Door Alarm | | ✓ | | | | |
| Lamp (Fluorescent) | | ✓ | | | | |

Ascent™ Storage Cabinet

Filtered Storage Solution

Engineering Drawing



1. Sample Port
2. Sentinel™ Control System

3. Electrical panel
4. Front panel


5. Secondary Filter
(For B, C, and E series only)

6. Blower
7. Main Filter

General Specifications, Ascent™ Storage Cabinet (4ft. series)

| Models | 220-240 VAC, 50-60 Hz, 1Ø | ASC-4A8 2040317 | ASC-4B8 2040284 | ASC-4C8 2040285 | ASC-4D8 2040325 | ASC-4E8 2040286 |
|---------------------------------|--|---|-----------------------------|-----------------------------|--------------------|-----------------------------|
| | 110-120 VAC, 50-60 Hz, 1Ø | ASC-4A9 2040318 | ASC-4B9 2040288 | ASC-4C9 2040289 | ASC-4D9 2040326 | ASC-4E9 2040290 |
| Nominal Size | 1.2 m (4') | | | | | |
| External Dimensions (W x D x H) | 1370 mm x 700 mm x 2050 mm (53.9" x 27.6" x 80.7") | | | | | |
| Door | Powder Coated electro-galvanized steel door with 6 mm thick Acrylic window and lockable Swing Handle | | | | | |
| Main Body | 1.2 mm 18 gauge electro-galvanized steel with white oven-baked epoxy-polyester ISOCIDE™ antimicrobial powder coated finish | | | | | |
| Number of Blower/s | 1 | | | | | |
| Air volume | 140 cfm | | | | | |
| Trays | Number of Trays | Standard: 8 trays (4 per section) Maximum: 12 trays (6 per section) | | | | |
| | Tray capacity per bottle size | For 2.5 -4L bottles: each tray can fit 9 bottles For 2.0 L bottles: each Tray can fit 12 bottles | | | | |
| | Maximum weight per tray | 40 kg | | | | |
| | Tray dimension | 600 mm x 600 mm (23.6" x 23.6") | | | | |
| | Tray material | Polypropylene | | | | |
| Filtration Elements | Main Filter | Nanocarb™ Filter (type A-H) | Nanocarb™ Filter (type A-H) | Nanocarb™ Filter (type A-H) | HEPA filter | HEPA filter |
| | Secondary Filter | None | HEPA filter | Nanocarb™ Filter (type A-H) | None | Nanocarb™ Filter (type A-H) |
| Control System | Sentinel™ Silver microprocessor control system | | | | | |
| VOC Sensor | ✓ | | | | | |
| Door Alarm | ✓ | | | | | |
| Lamp (Fluorescent) | ✓ | | | | | |

Engineering Drawing



- 1. Sample Port
- 2. Sentinel™ Control System
- 3. Electrical panel
- 4. Front panel
- 5. Secondary Filter
(For B, C, and E series only)
- 6. Blower
- 7. Main Filter

General Specifications, Ascent™ Storage Cabinet (6ft. series)

| | | | | | | |
|---------------------------------|--|---|-----------------------------|-----------------------------|--------------------|-----------------------------|
| Models | 220-240 VAC, 50-60 Hz, 1Ø | ASC-6A8 2040319 | ASC-6B8 2040321 | ASC-6C8 2040323 | ASC-6D8 2040327 | ASC-6E8 2040331 |
| | 110-120 VAC, 50-60 Hz, 1Ø | ASC-6A9 2040320 | ASC-6B9 2040322 | ASC-6C9 2040324 | ASC-6D9 2040328 | ASC-6E9 2040332 |
| Nominal Size | 1.8 m (6') | | | | | |
| External Dimensions (W x D x H) | 2055 mm x 700 mm x 2050 mm (80.9" x 27.6" x 80.7") | | | | | |
| Door | Powder Coated electro-galvanized steel door with 6 mm thick Acrylic window and lockable Swing Handle | | | | | |
| Main Body | 1.2 mm 18 gauge electro-galvanized steel with white oven-baked epoxy-polyester ISOCIDE™ antimicrobial powder coated finish | | | | | |
| Number of Blower/s | 2 | | | | | |
| Air volume | 140 cfm | | | | | |
| Trays | Number of Trays | Standard: 12 trays (4 per section) Maximum: 18 trays (6 per section) | | | | |
| | Tray capacity per bottle size | For 2.5 -4L bottles: each tray can fit 9 bottles For 2.0 L bottles: each Tray can fit 12 bottles | | | | |
| | Maximum weight per tray | 40 kg | | | | |
| | Tray dimension | 600 mm x 600 mm (23.6" x 23.6") | | | | |
| | Tray material | Polypropylene | | | | |
| Filtration Elements | Main Filter | Nanocarb™ Filter (type A-H) | Nanocarb™ Filter (type A-H) | Nanocarb™ Filter (type A-H) | HEPA filter | HEPA filter |
| | Secondary Filter | None | HEPA filter | Nanocarb™ Filter (type A-H) | None | Nanocarb™ Filter (type A-H) |
| Control System | Sentinel™ Silver microprocessor control system | | | | | |
| VOC Sensor | ✓ | | | | | |
| Door Alarm | ✓ | | | | | |
| Lamp (Fluorescent) | ✓ | | | | | |

Ascent™ Storage Cabinet

Filtered Storage Solution

CARBON FILTER OPTIONS

Dimension : 425 X 520 X 65 mm

| ITEM CODE | MODEL CODE | DESCRIPTION |
|-----------|------------|-----------------------------|
| 5080003 | CFA-11 | Standard filter |
| 5080095 | CFB-11 | Acid Filter |
| 5080028 | CFC-11 | Mercury Filter |
| 5080039 | CFD-11 | Sulphur Filter |
| 5080050 | CFE-11 | Halogen Filter |
| 5080061 | CFF-11 | Aldehyde Filter |
| 5080072 | CFG-11 | Ammonia / Amine Filter |
| 5080104 | CFH-11 | Ethers / Chloroforms Filter |

NUMBER OF CARBON FILTERS TO ORDER:

| Model | Qty (pcs) |
|--------------|-----------|
| ASC-2A/2B/2E | 1 |
| ASC-2C | 2 |
| ASC-4A/4B/4E | 1 |
| ASC-4C | 2 |
| ASC-6A/6B/6E | 2 |
| ASC-6C | 4 |

ACCESSORIES AND OPTIONS

| ITEM CODE | DESCRIPTION |
|-----------|---|
| 5021370 | Additional Polypropylene Tray 1 Set = 2 Trays: For 2 ft model – 1 set For 4 ft model – 2 sets For 6 ft model – 3 sets |

WE PROVIDE A ONE-STOP SOLUTION FROM INSTALLATION TO DECOMMISSIONING OF YOUR EQUIPMENT!

Esco Services takes professional pride in ensuring compliance with legal guidelines and ISO standards while helping your company minimize disruption or downtime of equipment; thus improving your company's productivity.



ACCREDITATION

Our Field Service Representatives have extensive annual up-to-date safety training. Our list of credentials includes NSF (for Biosafety cabinets) TÜV NORD (for microbiological & cytotoxic cabinets), NEBB (for clean room), IFBA (for biorisk management), and CETA-CNBT (for fume hoods).



STRONG GLOBAL NETWORK

We offer comprehensive and reliable services through our extensive global network of Field Service Representative to our clients that include a dependable inventory of ready parts, technical support, field-based engineering, factory engineers, and after-sales representatives.



FAST RESPONSE

On-site response within 24 hours for priority cases/ clients. Easy to reach customer service, through email and phone call.



GOOD GLOBAL REPUTATION

Our Service Team is dedicated to providing service of the highest standards. We provide timely, responsive service while meeting or even exceeding customer's expectations.

OUR WIDE RANGE OF QUALITY SERVICES



Installation and Operational Qualification

Esco Service Team provides installation services and ensures that the unit is operational upon installation.



Annual Recertification Testing and Calibration

Esco Field Service Representatives credentials include NSF (for Biosafety Cabinets), TÜV NORD (for Microbiological & Cytotoxic Cabinets), NEBB (for Clean Room), IFBA (for Biorisk Management), and CETA-CNBT (for Fume Hoods). Rest assured that your equipment will be properly certified according to international standards.



Decontamination

Esco offers safe and reliable decontamination services to ensure that your equipment is ready and safe to use.



Contract Manufacturing and Development Organization (CDMO)

Esco Aster, as a full-range services partner, supports its global clients from early stage discovery to market approval of clients' products. We provide process development (PD) and contract manufacturing (from upstream to downstream, formulation and filling), consultation, clinical research services bridging the translation gap from bench to bedside.



Repair, Maintenance and Replacement

Esco Service Engineers are the best in providing service for Esco products and other brands. Esco provides authentic parts for replacement and any components that might fail while the product is in service.

OTHER SERVICES OFFERED

With Esco service, customers can expect that breakdown calls are attended promptly and technical issues can be resolved faster than any other service provider.



Cleanroom Certification

Esco cleanroom services experts will verify your manufacturing environment ensuring performance requirement and the compliance of related industry standards.



BSL Validation

Esco is qualified to validate your laboratory according to your appropriate biosafety level. Facility design and operational parameters must meet the predetermined requirements by the WHO Biosafety manual prior to operation.



SCAN THIS QR TO KNOW MORE ABOUT OUR WORLD-CLASS SERVICE SOLUTIONS

ESCO

LIFESCIENCES GROUP

Improving Lives through Science

ESCO SCIENTIFIC



- Animal Research Workstation
- Biological Safety Cabinet
- CO₂ Incubator
- Ducted Fume Hood
- Ductless Fume Hood
- Filtered Storage Cabinet
- Laboratory Centrifuge
- Laboratory Oven and Incubator
- Laboratory Refrigerator and Freezer
- Laboratory Shaker
- Laminar Flow Cabinet
- PCR Cabinet
- Powder Weighing Balance Enclosure
- Thermal Cycler
- Ultra-low Temperature Freezer

ESCO MEDICAL



- Time-Lapse Incubator
- Benchtop Incubator
- ART Workstation
- CO₂ Incubator
- Anti-Vibration Table
- Gas Analyser

ESCO Healthcare

"Discovery to Delivery"



ESCO
PHARMA

ESCO
TaPestle

ESCO
VACCIXCELL


Pharmaceutical Manufacturing, Pharmacy Compounding,
and Bioprocessing Tools

ESCO ASTER



CRDMO Services

ESCO LIFESCIENCES GROUP
42 LOCATIONS IN 21 COUNTRIES ALL OVER THE WORLD



-  Global Offices
-  Licensee
-  Distributors
-  Factories
-  R&D Centers
-  Regional Distribution Centers

*Follow us on social media, download our apps,
and scan the QR code for more info.*



@EscoLifesciences



@EscoLifesciences



@EscoLifesci



@Esco



@EscoLifesciences



@EscoLifesciences



Esco Lifesciences



Esco Lifesciences



Esco Micro Pte. Ltd. • 19 Changi South Street 1 • Singapore 486779
Tel: +65 6542 0833 • mail@escolifesciences.com
www.escolifesciences.com

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA
Tel: +1 215-441-9661 • eti.admin@escolifesciences.com



Esco Lifesciences Group Offices: Bangladesh | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam