

Antimicrobial effectiveness evaluation of Isocide™ powder coating versus stainless steel plate

By: Dian Susanti and Alexander Atmadi

Abstract

Undesirable bioburdens can cause contamination on equipment and facilities used by laboratory professionals, where it is critical for them to keep everything clean and free from any health hazard. The bioburdens are present on many surfaces and it is hard to control their numbers effectively.

Some laboratory and cleanroom products are equipped with organism growth control mechanisms. One such mechanism is Isocide™ powder coating.

Two tests, based on JIS and ASTM standards, were performed to evaluate antimicrobial effectiveness of this powder coating. Stainless steel plate was used for this test as comparison. The test results showed that bacteria, yeast, and mold that came into contact with Isocide™ coated surfaces were inhibited and effectively eliminated.

Introduction

Over the years, numerous companies have been manufacturing controlled environment, laboratory and cleanroom equipment including biological safety cabinets, laminar airflow clean benches, animal handling workstations, laboratory fume hoods, ovens, incubators, PCR and IVF workstations. These products are widely used in many applications and facilities which are exposed to germs and microorganisms. Thus, there is a valid concern that contact with such objects may inflict illness from microorganisms such as bacteria, viruses, fungi.

It is very important to prevent biological contamination on the surfaces of laboratory equipment. Certain efforts have been undertaken to produce laboratory equipment with the ability to kill or inhibit the growth or reproduction of microorganisms – including the introduction of antimicrobial paint or powder coating.

There is no standard method published by either the EPA or BPD (EU Biocidal Products Directive) to determine the efficacy of antimicrobial paint. Many industry groups, such as ASTM, ISO, JIS, etc., publish their own standard methods that are primarily designed to determine the activity of antimicrobial agents on non-porous materials.

Below, two testing procedures were developed based on basic principles of JIS and ASTM testing standards, to provide both qualitative and quantitative data on the antimicrobial activity of Isocide™

Materials and methods

Microorganism Preparation.

Table 1. The tests used three types of microorganisms.

Table 1		
Bacteria	Yeast	Mold
<i>Bacillus subtilis var globigii</i>	<i>Saccharomyces cerevisiae</i>	<i>Aspergillus niger</i>
<i>Streptococcus epidermidis</i>	<i>Candida albicans</i>	
<i>Enterococcus faecalis</i>	<i>Rhodotula rubra</i>	
<i>Escherichia coli</i>		
<i>Staphylococcus aureus</i>		
<i>Serratia marcescens</i>		

Bacteria and yeast cultures were prepared by passing 2 activation phases, and their vegetative cells were prepared in phosphate buffer solution at pH 7.3 with 24 hours of growth for bacteria, and 36 hours for yeast cells. *Aspergillus niger* spore suspension was also used in this experiment with concentration 5 x 10⁸ spores/ml.

Testing procedures

A. Surface contact testing













1. Two type plates were used. Stainless steel and powder coating plates. Test plates were prepared from 50 x 50 mm electro galvanized steel plates, polished with zinc for rust protection, then sprayed with Isocide™ powder coating and baked in an oven at a temperature exceeding 180°C. These plates were swabbed with 70% Isopropyl alcohol to remove any visible dirt, left to dry and stored in the dark at 20°C prior to use.
2. All test microorganisms were pipetted with 0.1 ml pipettes into agar plates, and spread to the entire agar surface. Test plates were put in the center of inoculated agar. Bacteria cultures were spread into Trypticase soy agar, yeast culture inoculated to Saboraud dextrose agar, and mold was spread to Potato dextrose agar surface.
3. Growth of test microorganisms was observed 24-48 hours for bacteria, 48-72 hours for yeast, while mold was observed until 1 week of incubation period. Any growth between agar and test plate surface contact was documented.

B. Determination of log reduction of bioburden

1. Test plates were prepared from 20 x 50 mm electro galvanized steel plates, polished with zinc rust protector, and coated with Isocide™ powder coating and baked at a temperature exceeding 180°C. Those plates were swabbed with 70% Isopropyl alcohol to remove any visible dirt, left to dry and stored in the dark at 20°C.
2. A stainless pan was prepared to hold the test plates during contact time, and this pan was equipped with organic decontaminator sheet, to remove any contaminant at the back surface of test plates.
3. Test plates were arranged accordingly on top of the decontaminator sheet in the pan, and pipetted with 0.1 ml of test microorganisms, and then the pan was covered with plastic wrap to keep the plates safe from unnecessary materials and contaminants. These plates were put inside incubator without humidity control at 25°C.
4. Tested plates were sampled at 24 and 48 hours of contact time, the remaining CFU (colony forming unit) of test microorganisms were quantified by serial dilution method, and log reduction was retrieved by comparing initial concentration with CFU remains in the test plates.

Results

The resulting different growth responses between the colonies with direct contact (under the test plate), minimum contact (on the plate's edge), and no contact (outside the plate) are shown in Figure 1 and described in Table 2.

Agar control	Agar with plate	Agar with SS plate
		
<i>Bacillus subtilis</i> var <i>globigii</i>		
		
<i>Escherichia coli</i>		
		
<i>Enterococcus faecalis</i>		
		




Figure 1. Test plates from Test Procedure

Table 2: Results from Test Procedure A

Microorganisms	Growth observation		
	Isocide plate	SS plate	Outside of plate
<i>Bacillus subtilis</i> varglobigii	+	++	++++
<i>Serratiamarcescens</i>	-	++	++++
<i>Staphylococcus aureus</i>	-	+	++++
<i>Streptococcus epidermidis</i>	-	++	++++
<i>Enterococcus faecalis</i>	-	++	++++
<i>Escherichia coli</i>	-	++	++++
<i>Saccharomyces cerevisiae</i>	-	++	++++
<i>Candida albicans</i>	-	++	++++
<i>Rhodotularubra</i>	-	++	++++
<i>Aspergillusniger</i>	-	+++	++++

(-) No growth; (+ to +++) Representing more or less dense population of microorganisms

Table 3: Results from Test Procedure B - Microbial log reduction after incubation with Isocide™ powder coating. Log reduction is calculated as the difference between log initial concentration and log final concentration

Microorganisms	Initial concentration (CFU/ml)	Cells' logarithmic value	24 hr		48hr	
			CFU remain	Log reduction	CFU remain	Log reduction
<i>Bacillus subtilis</i> varglobigii	3.0×10^5	5.48	57	3.70	36	3.90
<i>Serratiamarcescens</i>	2.4×10^5	5.38	0	5.38	0	5.38
<i>Staphylococcus aureus</i>	3.3×10^5	5.52	45	3.87	17	4.29
<i>Streptococcus epidermidis</i>	1.2×10^5	5.08	40	3.48	15	3.90
<i>Enterococcus faecalis</i>	1.6×10^5	5.20	0	5.20	0	5.20
<i>Escherichia coli</i>	1.3×10^5	5.11	5	4.50	0	5.11
<i>Saccharomyces cerevisiae</i>	9.0×10^4	4.95	10	3.95	0	4.95
<i>Candida albicans</i>	1.0×10^4	4.00	50	2.30	12	2.90
<i>Rhodotularubra</i>	4.0×10^4	4.60	90	2.65	30	3.12




Figure 2. Isocide™ performance between different microorganisms

Observations

1. Table 2 show no microorganism growth, except for *B. subtilis*, in areas where Isocide™ was present.
2. During the time of observation, some *B.subtilis* colonies were growing under the test plate as described at figure 1.
3. Isocide™ performance varies, depending upon the type of microorganism. The best performance showed on
4. *Serratia marcescens* and *Enterococcus faecalis* with no single colony growth after one day incubation
5. Table 2, show that all test microorganism capable growth under the SS plates. This indicating that SS cannot act as anti-microbial as powder coating plates.

Conclusions

This study supports the idea that the addition of antimicrobial coating can be effective in preventing or reducing microbial growth. The antimicrobial evaluations performed in this study were based on the zone of inhibition and JIS 2801 tests. We found that Isocide™ powder coating perfectly can serve as an antimicrobial agent compared with stainless steel plates since it does effect the microbial population growth.

References

- [1] Barret, L., Emerentiana, S., Darrel, S. 2007. Antimicrobial Coating. <http://www.battelle.org/forecast/home.stm>.
- [2] Sadasivan, L. 2007. Antimicrobial coating. <http://www.battelle.org/forecast/driver.stm>. Siva Microbiological Solutions LLC, Bristol, PA
- [3] Tesh, E.M. The Regulation of Antimicrobials in Paints and Surface Treatments. Paint and Coatings Industry 2004, 7