



PIPETTE A2,
L-2107

PIPETTE A10,
L-2108

PIPETTE A20,
L-2109

PIPETTE
A100, L-2111

PIPETTE
A200, L-2112

PIPETTE A1000,
L-2113

PIPETTE
A5000,
L-2114

PIPETTE
A10ML,
L-2115

Mechanical Pipette

For precise measurement and productive liquid sample transfer



MECHANICAL PIPETTES



Ushering in a new age of traceability research, the A-Series Mechanical Pipette is a new line of pipettes aimed at improving ergonomics and comfort while maintaining its accuracy and reliability.

Key Features

Control Button

- Very low operating force
- Colour-coded based on volume
- Positioned for perfect ergonomics

Volume Adjustment

- Can turn to adjust volume

Ejector

- Low tip ejection force of less than or equal to 4 newtons
- Positioned for perfect ergonomics

Volume Display


- 4-digit display
- Magnified for comfortable viewing

Quick Connection Tip

- Remove lower part quickly

Tip Cone

- Spring-loaded tip cone for ease of connecting and ejecting pipette tips



1. Control Button
2. Finger Rest

3. Quick Connection Tip
4. Volume Adjustment

5. Ejector
6. Volume Display

7. Tip Cone

Product Features

- 4-digit volume display
- Magnified display for comfortable viewing
- Chemical resistant
- Up to 121°C sterilization
- Very low operating force
- Colour-coded according to volume
- Positioned for perfect ergonomics
- Turn to adjust volume
- Antimicrobial surface
- No need for greasing of the pipette piston after autoclaving

Mechanical Pipette Options and Accessories

Item Code	Model Code	Description	Accuracy	Repeatability
5171116	PIPETTE A2, L-2107	Volume Range of 0.1 - 2 µl, 1 Pack (20 pcs per pack)	±12.0-1.5%	≤6.00-0.70%
5171117	PIPETTE A10, L-2108	Volume Range of 0.5 - 10 µl, 1 Pack (20 pcs per pack)	±2.5-1.0%	≤1.50-0.40%
5171118	PIPETTE A20, L-2109	Volume Range of 2 - 20 µl, 1 Pack (20 pcs per pack)	±2.5-1.0%	≤1.50-0.30%
5171119	PIPETTE A100, L-2111	Volume Range of 10 - 100 µl, 1 Pack (20 pcs per pack)	±1.8-0.8%	≤0.50-0.15%
5171120	PIPETTE A200, L-2112	Volume Range of 20 - 200 µl, 1 Pack (20 pcs per pack)	±1.8-0.6%	≤0.50-0.15%
5171121	PIPETTE A1000, L-2113	Volume Range of 100 - 1000 µl, 1 Pack (20 pcs per pack)	±1.5-0.6%	≤0.30-0.15%
5171122	PIPETTE A5000, L-2114	Volume Range of 500 - 5000 µl, 1 Pack (20 pcs per pack)	±1.2-0.6%	≤0.30-0.15%
5171123	PIPETTE A10ML, L-2115	Volume Range of 1 - 10 mL, 1 Pack (20 pcs per pack)	±2.5-0.6%	≤0.30-0.15%

Complete Set


Item Code	Description
5171125	Complete Set of 8 Pipettes (A2, A10, A20, A100, A200, A1000, A5000, A10ML) & 1 Pipette Stand (HW-9)

Accessories

Item Code	Model Code	Description
5171124	PIPETTE STAND HW-9, L-2039A	Pipette stand that holds up to 9 Pipettes, 1 Pack (20 pcs per pack)

ESCO LIFESCIENCES GROUP

42 LOCATIONS IN 21 COUNTRIES ALL OVER THE WORLD



- Global Offices
- Distributors
- Factories
- R&D Centers
- Regional Distribution Centers

*Follow us on social media, download our apps,
and scan the QR code for more info.*



ESCO[®]

LIFESCIENCES GROUP

Esco Micro Pte. Ltd. • 19 Changi South Street 1 • Singapore 486779
Tel +65 6542 0833 • mail@escolifesciences.com
www.escolifesciences.com

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA
Tel: +1 215-441-9661 • eti.admin@escolifesciences.com

Esco Lifesciences Group Offices: Bangladesh | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania
| Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam