




## **Automatic Temperature Measurement and Disinfection Machine**



## Automatic Temperature Measurement and Disinfection Machine

A contactless multifunctional device that can take quick and accurate temperature measurements and simultaneously dispense liquid sanitizer. It is equipped with an automatic alarm that triggers when the measured body temperature exceeds the preset alarm temperature.



### Large Capacity

Comes with 1,000 mL liquid soap capacity.



### Automatic Induction

Convenient to use with intelligent point saving and automatic wake-up function.



### Temperature Measurement

Can take quick temperature measurements.



### Alarm Function

Equipped with visual and audible alarm to signal high temperature.



### °C/°F Unit Available

Can display temperature in °C or °F.



### High Accuracy

High-precision sensor with stable and reliable performance.

Align the palm of your hand with the temperature reading sensor and liquid outlet to allow the device to measure the temperature and automatically dispense the liquid sanitizer.



User Guide

USB cable & charger

Screws & rubber plugs









#### Automatic Temperature Measurement and Disinfection Machine

Accuracy	±0.2°C (34 - 45°C, put it in the operating environment for 30 minutes before use)
Alarm type	Beeping and flashing LED indicator
Measuring distance	5 - 10 cm
LED display	Digital
Interface	USB Type-C Port
Bottle capacity	1000 mL
Mounting type	Wall-hanging, Tripod-fixed
Working temperature	8.10 - 40°C (15 - 35°C Recommended)
Measuring range	0 - 50°C
Standby time	5 seconds for non-working, dynamic loop with white dot
Dimensions (W x D x H)	119 x 133 x 280 mm
Net weight	702 g
Input	DC 5V 1A
Power consumption	3W
Liquid type	Hand sanitizer


# ESCO LIFESCIENCES GROUP

42 LOCATIONS IN 21 COUNTRIES ALL OVER THE WORLD



-  Global Offices
-  Joint Ventures
-  Licensee
-  Distributors
-  Factories
-  R&D Centers
-  Regional Distribution Centers

Join us on Social Media and Download our Apps!



Filtra Check App  
The World's First Chemical  
Evaluation App



Esco Social Media  
and Apps

# ESCO

LIFESCIENCES GROUP

Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777  
Tel: +65 6542 0833 • Fax: +65 6542 6920 • mail@escolifesciences.com  
www.escolifesciences.com

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA  
Tel: +1 215-441-9661 • Fax 484-698-7757  
eti.admin@escolifesciences.com

Esco Lifesciences Group Offices: Bangladesh | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania  
| Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam