


MARK	REVISION	DATE	REV BY	ACC
A	FIRST RELEASE DRAWING	18/10/2022	Firman	Indra S.



GENERAL SPECIFICATION	
Model	LA2-SS9 G4 10" with STA-5A0
External Dimensions (WxDxH)	1645x873x2285-2485 mm(64.7"x34.4"x90.0"-97.8")
External Dimensions with STA-5A0 (WxDxH)	1679x873x2285-2485 mm(66.1"x34.4"x90.0"-97.8")
Gross Internal Dimensions (WxDxH)	1525x825x720 mm(60.0"x24.6"x28.3")
Usable Work Area	0.80 m² (8.5 sq.ft)
Working Opening	254 mm (10")
Maximum Opening	575 mm (22.6")
Filter Typical Efficiency	>99.999% for particle size between 0.1 to 0.3 microns, per IEST-RP-CC001.3
Light Intensity (Dimmable)	>1200 Lux (>111 foot-candles)
Cabinet Construction	Main Body: 1.5 mm(0.06")16 gauge EG Steel With Isocide Oven-Baked Epoxy-Polyester Powder Coating Work Zone: Stainless Steel 304, 48 Finish. Tray: 1.5mm(0.06")16 gauge. Walls: 1.5 mm(0.06")16 gauge
Electrical	Power Rating: 9
Average Airflow	0.53 ± 0.025 m/s (105 fpm ± 5 fpm)
Velocity	Downflow 0.33 ± 0.025 m/s (65 fpm ± 5 fpm)

Notes (PLEASE READ!) -All drawings and specifications must be signed and approved by the end customer. -The factory will not proceed with design and fabrication without the signature of approval. -Contact ESCO immediately if there are any changes or clarifications needed. -Design changes may result in pricing changes as well. -Any requests of revisions after the final approval will not be entertained.		NOTE		THIRD ANGLE VIEW		DIMENSION: MM / INCH		SCALE: NTS	
		End Customer :							