



Leading Global Pharmaceutical Company Commissions Esco GB Ltd for Turbulent Flow Containment Isolator

An Esco **Containment Barrier Isolator – Turbulent Flow (CBI-T)** will be installed in the client's facility in Scotland. The CBI-T model was chosen as it provides the perfect solution to meet the client's process requirements.

The CBI-T is required as their process is to accurately weigh and dispense small quantities of Active Pharmaceutical Ingredients (APIs) during quality control tests. The rear wall of the unit is equipped with a screen for weight display and other operational functions, making the CBI-T well suited to this process. It is also integrated with controllable lighting for light-sensitive products.

The CBI-T is designed to facilitate the isolation of a potent product or process while providing operator, external environmental, and product protection.

For a more ergonomic solution for various operators who will use the equipment, the unit is equipped with an electric height-adjustable stand so that the operator can select a setting suited to their preference.

By early 2020, Esco Pharma will also deliver to the same client a customised downflow booth, equipped with material and personnel airlock for bulk weighing and dispensing operations.

For more information please visit www.escopharma.com.

Esco Pharma

21 Changi South Street 1 Singapore 486777 • Tel: +65 65420833 • Fax: +65 65426920

Email: csis.pharma@escoglobal.com • Website: www.escopharma.com

Esco Global Offices: Bangladesh | Cameroon | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam

ESCO

WORLD CLASS. WORLDWIDE.